

Gluteal Smooth

Medium, High Projection

Breast & Gluteal Implant
BellaGel






· *Gluteal implants are used for augmentation, contouring or symmetry correction after trauma and implant exchange. More cohesive & harder gel filled to withstand the patient's weight and movement.*

Specification

Picture	Shape	Surface	Height	Projection	Volume	Cat. No.	
						Smooth	Texture
	Round	Smooth	N/A	Medium	160~340cc	GRSZ-M cc	N/A
				High	220~430cc	GRSZ-H cc	

Gluteal Smooth

Medium, High Projection



Medium Projection

Cat. No	Volume (cc)	Diameter (cm)	ProjectionA (cm)	ProjectionB (cm)
GRSZ-M 160	160	10.0	3.0	3.4
GRSZ-M 230	230	11.0	3.4	3.8
GRSZ-M 265	265	12.0	3.5	3.9
GRSZ-M 340	340	13.0	3.8	4.2

High Projection

Cat. No	Volume (cc)	Diameter (cm)	ProjectionA (cm)	ProjectionB (cm)
GRSZ-H 220	220	10.0	4.1	4.4
GRSZ-H 280	280	11.0	4.3	4.6
GRSZ-H 340	340	12.0	4.5	4.9
GRSZ-H 430	430	13.0	4.8	5.2